

Test results

Corresponding to our wide experience in testing products we know that the continuous contact submerged in a 10% solution of the cleaning agent EW80 clean and the disinfectant EW80 des (instead of only 1 % as VAH / DGHM certified) for about 4 weeks simulates at least 100 cycles.

In addition we put the masks in a washing machine with 50° C to test the mechanical and thermal stress. In between and after the end of the test we could not assert any visual or tactile difference in comparison to the not tested items.

All these careful and systematic examinations prove, that the Laerdal masks are compatible to the EW80 products.

Dortmund, 15.02. 2016



.....
General Manager



.....
Test performer