



EW80 Systeme GmbH
Geschäftsführer Hanno von Roth
HRB 4009 AG Dortmund

Untere Brinkstraße 80
D-44141 Dortmund
Telefon (02 31) 94 10 14-0
Telefax (02 31) 94 10 14-10
www.ew80.de
dekopause.ew80.de

Test results

Corresponding to our wide experience in testing products we know that the continuous contact submerged in a 10% solution of the cleaning agent EW80 clean and the disinfectant EW80 des (instead of only 1 % as VAH / DGHM certified) for about 4 weeks simulates at least 100 cycles.

In addition we put the masks in a washing machine with 10% solution and for 20 cycles by 50° C to test the mechanical and thermal stress. The 10% solution instead of 1% of EW80 clean and EW80 des simulates another 100 cycles for each product. In between and after the end of the test we could not assert any visual or tactile difference in comparison to the not tested items.

All these careful and systematic examinations prove, that the Ambu masks are compatible to the EW80 products.

Dortmund, 23.02.2016

.....
General Manger

.....
Test performer